



How To Add In Logbook


1 Click "Setup Page"



2 Click "Location Setup"




3 Click "Add Location"



The screenshot shows the 'guard' application interface. At the top, there is a navigation bar with a user profile icon, the ID 'IG02362', 'Contact Support', 'Welcome Donald', a notification bell, a settings gear, and a 'Log Out' link. Below the navigation bar, the page title is 'Setup' and the breadcrumb is 'Location Setup'. A search bar labeled 'Search: Search Locations' is present. The main content area is a table with a header 'Actions' and two columns of blue icons: a pencil for edit and a trash can for delete. The 'Add Location' button in the top right corner is circled in orange.

4 Click and enter the location name in the box.




The screenshot shows the 'Location Setup' page with an 'Add New Location' modal dialog open. The dialog has a title 'Add New Location' and a label 'Location Name:' followed by a text input field containing the placeholder text 'Enter Location'. The input field is circled in orange. Below the input field are 'Close' and 'Save' buttons. The background page is dimmed, showing the same table and search bar as in the previous screenshot. At the bottom of the page, there is a footer with 'Items per page: 10' and '1 - 10 of 54' along with navigation arrows.

5


Click "Save"

The screenshot shows a web application interface for 'Location Setup'. At the top, there is a search bar labeled 'Search: Search Locations'. Below it is a table with columns for location names and actions. The table contains several rows, with the first row partially visible as 'partment A'. To the right of the table, under the heading 'Actions', there are two columns of icons: a pencil icon for editing and a trash can icon for deleting. A modal dialog titled 'Add New Location' is open in the center. It has a 'Location Name:' label and a text input field containing '1 north'. At the bottom of the modal, there are two buttons: 'Close' and 'Save'. The 'Save' button is highlighted with an orange circle. At the bottom of the page, there is a footer with the text 'elliguard'.



How to Edit Locations

- 1 Search for a location in the search field.




The screenshot shows the Intelliguard web interface. At the top, there is a navigation bar with the Intelliguard logo, user ID IG02362, Contact Support, Welcome Donald, and Log Out. Below the navigation bar is the "Location Setup" page header with an "Add Location" button. A search bar labeled "Search: Search Locations" is highlighted with an orange circle. Below the search bar is a table with the following data:

Location Name	Actions
1 north	 
3M	 
5 North	 
6 North	 
6 West	 
Cardiac Cath Lab	 
CATH LAB	 
Central Supply	 
CT 1 (new)	 
CT 2 (old)	 

At the bottom of the table, there is a pagination control showing "Items per page: 10" and "1 - 10 of 55".


- 2 Click the "Edit" button under actions



The screenshot shows the Intelliguard web interface, identical to the first screenshot. In this view, the "Edit" button (pencil icon) under the "Actions" column for the first row, "1 north", is highlighted with an orange circle.


3 Enter the new name and click "Save"

The screenshot displays a web application interface for "Location Setup". At the top, there is a search bar labeled "Search: Search Locations". Below this is a table with a header "Actions" containing edit and delete icons for each row. A modal dialog titled "Edit Location Name" is open in the center. It features a text input field labeled "Location Name:" with the text "1 north west" entered. At the bottom of the modal, there are two buttons: "Close" and "Save". The "Save" button is highlighted with an orange circle. At the bottom of the page, there is a footer with the text "elliguard".



How to Delete a Location

1 Click the "trash can" under actions



2 It will ask you to confirm that you want to delete the location and click "Yes"

