



How to Create a New Cart Type

1 Click "Setup Page"



2 Click "Cart Setup"



3

Click "Add New Cart Type"

The screenshot shows the 'guard' application interface. At the top, there is a navigation bar with the logo 'guard', a user profile 'IG02362', 'Contact Support', 'Welcome Donald', and 'Log Out'. Below the navigation bar, there is a search bar labeled 'Search: Search Cart Types'. To the right of the search bar, there is a blue button labeled 'Add New Cart Type' which is circled in orange. Below the search bar, there is a table with columns for 'Actions' and several rows of data. Each row contains three blue icons: a pencil, a camera, and a hash symbol.

4

Enter the new Cart Type name and then click "Save"

The screenshot shows the 'guard' application interface with a modal dialog box open. The dialog box is titled 'Add New Cart Type' and contains a text input field with the value '1 new cart'. Below the input field, there are two buttons: 'Close' and 'Save'. The 'Save' button is circled in orange. The background of the application is dimmed, showing the same search bar and table as in the previous screenshot. At the bottom of the application, there is a footer with the text 'Items per page: 10' and '1 - 10 of 13'.

How to Add Cart Numbers

- 1 Click the "#" under actions next to the correct cart name.

Intelliguard IG02362 Contact Support Welcome Donald Log Out

Cart Setup

Logbook / Setup / Cart Setup

Search: Search Cart Types

Cart Type ↑	Actions
1 new cart	
Adult Crash Cart	
Ancillary Supply Cart	
Ben's Cart	
Diaper Cart	
Example Cart	
Helicopter Cart	
My New Cart	
Neonatal Crash Cart	
Nutraceutical Cart	

Items per page: 10 1 - 10 of 14

© Copyright 2025 Intelliguard

- 2 Click in the box and enter in the cart numbers, then click "save"

Intelliguard IG02362 Contact Support Welcome Donald Log Out

Cart Setup

Logbook / Setup / Cart Setup

Search: Search Cart Types

Manage Cart Numbers - 1 new cart

Click on a panel below to **Add** or **Edit** cart numbers for **Rady Children's Hospital**

Add New Cart Numbers Below:

New Cart Number: Save

Edit Existing Cart Number(s) Below:

Edit Cart Numbers associated with this Cart Type: Save

Cancel Done

Cart Type ↑	Actions
1 new cart	
Adult Crash Cart	
Ancillary Supply Cart	
Ben's Cart	
Diaper Cart	
Example Cart	
Helicopter Cart	
My New Cart	
Neonatal Crash Cart	
Nutraceutical Cart	

How to Edit Cart Numbers

1 Click the "#" under actions next to the correct cart nam

Intelliguard IG02362 Contact Support Welcome Donald Log Out

Cart Setup

Logbook / Setup / Cart Setup

Search: Search Cart Types

Cart Type ↑	Actions		
1 new cart			
Adult Crash Cart			
Ancillary Supply Cart			
Ben's Cart			
Diaper Cart			
Example Cart			
Helicopter Cart			
My New Cart			
Neonatal Crash Cart			
Nutraceutical Cart			

Items per page: 10 1 - 10 of 14

© Copyright 2025 Intelliguard

2 Click into the edit box.

Intelliguard IG02362 Contact Support Welcome Donald Log Out

Manage Cart Numbers - 1 new cart

Click on a panel below to **Add** or **Edit** cart numbers for **Rady Children's Hospital**

Add New Cart Numbers Below:

Existing Cart Numbers:
001, 002, 003, 004, 005, 007, 008

New Cart Number: **Save**

Edit Existing Cart Number(s) Below:


Edit Cart Numbers associated with this Cart Type:

001	002	003
004	005	007
008		


Save

Cancel **Done**

3 Click the number you want to edit and change it.




4 Click "Save"




How to Add Trays to a Cart

1 Click the briefcase button under actions.



2 Click the checkbox next to the appropriate trays that you want to add.



3 Click "Save"

Associate Tray Types with: 1 new cart


Check the following trays that you would like associated with this Cart Type:

AAA Tray	<input checked="" type="checkbox"/>
AC Tray	<input type="checkbox"/>
Adrian's Tray	<input checked="" type="checkbox"/>
Adult Anesthesia and Pain Tray	<input type="checkbox"/>
Adult Code Cart Tray	<input type="checkbox"/>
Adult Crash Cart Tray	<input type="checkbox"/>
Delete Test	<input checked="" type="checkbox"/>
Demco Anesthesia Tray	<input type="checkbox"/>

Close Save


Items per page: 10 1 - 10 of 14

elliguard




How to Edit Cart Name

- 1 Click the "Edit" button under the actions.



The screenshot shows the Intelliguard 'Cart Setup' interface. At the top, there is a navigation bar with the Intelliguard logo, user ID 'IG02362', 'Contact Support', and a 'Welcome' message. A green notification banner at the top right says 'Cart Type Updated.' Below the navigation bar, the page title 'Cart Setup' is displayed, along with a breadcrumb trail 'Logbook / Setup / Cart Setup' and an 'Add New Cart Type' button. A search bar labeled 'Search: Search Cart Types' is positioned above a table of cart types. The table has two columns: 'Cart Type' and 'Actions'. The 'Cart Type' column lists various cart types, including '1 new cart', 'Adult Crash Cart', 'Ancillary Supply Cart', 'Ben's Cart', 'Diaper Cart', 'Example Cart', 'Helicopter Cart', 'My New Cart', 'Neonatal Crash Cart', and 'Nutraceutical Cart'. The 'Actions' column contains three icons: a pencil (edit), a camera (clone), and a hash symbol (#). The 'Edit' icon for the '1 new cart' row is circled in orange. At the bottom of the table, there is a pagination control showing 'Items per page: 10' and '1 - 10 of 14' items, along with navigation arrows. A copyright notice '© Copyright 2025 Intelliguard' is visible at the bottom left.

- 2 Here you can change the name of the Cart and then click "Save".



This screenshot shows the same 'Cart Setup' page as above, but with a modal dialog box open. The dialog box is titled 'Edit Existing Cart Type: 1 new cart' and contains a text input field labeled 'Current Cart Type:' with the value '1 new carta'. Below the input field are two buttons: 'Close' and 'Save'. The 'Save' button is circled in orange. The background of the page is dimmed, and the 'Edit' icon from the previous screenshot is still visible. The pagination control at the bottom shows 'Items per page: 10' and '1 - 10 of 14' items.