



# How to Encode Pharmacy Inventory

1 Click "Encode Inventory"



2 Click this field, then Scan the barcode on the medication.



3

Type the Lot number into the field and then it will have you confirm it and type the lot number in again.

Inventory

Inventory

Confirmation

Product Code: 00143-9924-90  
Item Type: Medication  
Item: ceFAZolin 1 GM Inj (West-Ward Pharmaceuticals Corp.)  
Manufacturer: West-Ward Pharmaceuticals Corp.  
Package Type: Vial

Back

Lot Number: ZP521

Confirm Lot Number: |

Expiration Date: mm/dd/yyyy

Confirm Expiration Date: mm/dd/yyyy

Place items inside of enclosure, close and secure the door, then press scan.

Items Scanned: Scan

Undo and Save

Undo all fields and clear form. Clear

Acknowledge and clear form. Accept

4

Click the "mm/dd/yyyy" field and enter in the expiration date

Inventory

Inventory

Confirmation

Product Code: 00143-9924-90  
Item Type: Medication  
Item: ceFAZolin 1 GM Inj (West-Ward Pharmaceuticals Corp.)  
Manufacturer: West-Ward Pharmaceuticals Corp.  
Package Type: Vial

Back

Lot Number: ZP521

Confirm Lot Number: ZP521

Expiration Date: mm/dd/yyyy

Confirm Expiration Date: mm/dd/yyyy

Place items inside of enclosure, close and secure the door, then press scan.

Items Scanned: Scan

Undo and Save

Undo all fields and clear form. Clear

Acknowledge and clear form. Accept

5

It will have you confirm and type in the expiration date again

Confirmation

Product Code: 00143-9924-90  
Item Type: Medication  
Item: ceFAZolin 1 GM Inj (West-Ward Pharmaceuticals Corp.)  
Manufacturer: West-Ward Pharmaceuticals Corp.  
Package Type: Vial

Back

Lot Number: ZP521  
Confirm Lot Number: ZP521  
Expiration Date: 5/31/2025  
Confirm Expiration Date: mm/dd/yyyy

Place items inside of enclosure, close and secure the door, then press scan.

Items Scanned:  Scan

Undo all fields and clear form. Acknowledge and clear form.

Clear Accept


guard

6

Place all of the same medications with the same lot number inside the enclosure



Be sure you know how many medications are being put inside of the enclosure and that the Lot number and expiration date are the same on all of them.



## 7 Click "Scan"

Confirmation

Product Code: 00143-9924-90  
Item Type: Medication  
Item: ceFAZolin 1 GM Inj (West-Ward Pharmaceuticals Corp.)  
Manufacturer: West-Ward Pharmaceuticals Corp.  
Package Type: Vial

Back

Lot Number: XP523  
Confirm Lot Number: XP523  
Expiration Date: 5/31/2025

Place items inside of enclosure, close and secure the door, then press scan.

Items Scanned:  Scan

Undo all fields and clear form. Acknowledge and clear form.

Clear Accept

guard

## 8 The number of items scanned will populate, double check that this number matches what you originally counted.

## 9 Click "Accept"

Confirmation

Product Code: 00143-9924-90  
Item Type: Medication  
Item: ceFAZolin 1 GM Inj (West-Ward Pharmaceuticals Corp.)  
Manufacturer: West-Ward Pharmaceuticals Corp.  
Package Type: Vial

Back

Lot Number: XP523  
Confirm Lot Number: XP523  
Expiration Date: 5/31/2025

Place items inside of enclosure, close and secure the door, then press scan.

Items Scanned: 1 Scan

Undo all fields and clear form. Acknowledge and clear form.

Clear Accept

guard